

KOMPACT



Model Shown : 235-10 + accessories



Examination table
235-10 / 235-15

Product Identification

Model Number :

Serial Number :

Date of Purchase :

Name of Owner / Facility :

Name of Dealer :

Dealer's Phone Number :

Promotal Authorized Service Company :

Legal Notice


PROMOTAL
22, rue de Saint-Denis de Gastines
B.P. 26 - 53500 ERNÉE Cedex
FRANCE
Tél. : +33 (0)2 43 05 12 70
Fax : +33 (0)2 43 05 68 99
internet : www.promotal.com

The descriptions and specifications contained in this Operating Manual are deemed correct at the time of printing.

Promotal, however, reserves the right to modify its models and its procedures or render them obsolete without notice.

Before any order, we recommend that our customers consult a local sales manager.

TYPE PLATE



Year(YY)/Week(WW) of manufacture



Important information	6
<i>Safety symbols</i>	6
<i>Applied parts</i>	7
<i>Electromagnetic interference</i>	7
Unpacking precautions	8
<i>Medical device delivered on a wooden pallet</i>	8
<i>Storage conditions</i>	8
<i>Conditions of use</i>	8
<i>Unpacking and Installation</i>	8
<i>Step by step</i>	8
<i>Check</i>	9
<i>Remark:</i>	9
<i>Authorised EC Representative</i>	9
Cleaning protocol	10
<i>Warning</i>	10
<i>Cleaning/Disinfecting</i>	10
User manual – Kompact	11
<i>Intended purpose</i>	11
<i>Patient weight capacity</i>	11
<i>Protection against penetration of liquids</i>	11
<i>Characteristics</i>	11
<i>Dimensions / Installation precautions</i>	12
<i>235-10 / 235-15</i>	12
<i>235-10 / 235-15 + Option 272-01</i>	12
<i>Laying out the sleeping surface</i>	13
<i>One piece upholstery</i>	13
<i>Legrest upholstery</i>	14
Safety	15
Using the Kompact	16
<i>Backrest tilt adjustment</i>	16
<i>Trendelenburg position</i>	16
<i>Principle</i>	16
<i>Leg rest</i>	17
<i>Using the leg rest</i>	17
<i>Footstep</i>	17
<i>Using the footstep</i>	17
<i>Paper Roll</i>	18
<i>Paper roll installation</i>	18
<i>Using the paper roll</i>	18
<i>Gynaecology debris tray</i>	19
<i>Using the gynaecology debris tray</i>	19
<i>Integrated stirrups</i>	20
<i>Installing the integrated stirrups</i>	20

Summary

The levellers	20
<i>Adjustment of the levellers</i>	20
Drawers	21
Braking	22
<i>Using the brakes</i>	22
Accessories	23
Rear tray for storing paper rolls	23
<i>Set up</i>	23
2 removable plastic trays / Side drawers	23
<i>Set up</i>	23
Pair of stainless steel rails 25x10	24
<i>Remove plastic end caps</i>	24
<i>install the pair of rails</i>	24
Sliding clamp for 25 x 10 rail	25
<i>(reference 878-12)</i>	25
<i>Installing the clamp</i>	25
<i>Using the clamp</i>	25
Rotating clamp for 25 x 10 rail	26
<i>(reference 879-10)</i>	26
<i>Installing the clamps</i>	26
Auto-blocking I.V. pole	27
<i>(reference 2985-01)</i>	27
<i>Installing the IV pole</i>	27
<i>2 hook IV pole</i>	27
<i>(reference 985-01)</i>	27
<i>Installing the IV pole</i>	27
<i>Adjusting the height of the IV pole</i>	27
Arm Rest	28
<i>Adjusting the arm rest</i>	28
<i>Pair of legrests GOEPEL</i>	28
<i>Installing the legrests</i>	28
Lamp holder	29
<i>(reference 295-01)</i>	29
<i>Fixing of the lamp on the support</i>	29
Adjustable paper-roll holder	30
<i>(reference 2056-02)</i>	30
<i>Installating the paper roll holder</i>	30
<i>Installation of the paper roll</i>	30
Leg rest flush with seat	31
<i>(reference 272-01)</i>	31
<i>Using the leg rest</i>	31
Additional mains plug and equipotential terminal	32
<i>(reference 240-01EU / 240-01UK / 240-01US)</i>	32

<i>Using the additional mains plug</i>	32
<i>Equipotential terminal</i>	33
Notes	34
Lifespan of the Medical Device	35
<i>The wear parts are:</i>	35
Compulsory / specific maintenance	35
<i>Once a year, ensure that the following checks are carried out by a qualified technician (contact your dealer):</i>	35
<i>Once a month :</i>	35
Medical device end of service life	35
Fuse replacement	36
<i>(Only if option 240-01EU or 240-01UL or 240-01US is present)</i>	36
Maintenance notebook	37
Warranty Information	38
Warranty	38
<i>Obligations</i>	38
<i>Exclusions</i>	38
<i>Exclusive Remedy</i>	38
<i>No Authorization</i>	38

Important information

Safety symbols



Warning sign

Information marked with this symbol must be read and strictly complied with!

Remark

Draws your attention to a procedure, practice or situation.



Humidity limitation



Maximum number of stacked pallets



Correct orientation for transportation



Atmospheric pressure limitation



Temperature limitation



Fragile



Keep in a dry place



Not to use sharp objects



Manufacturer



CE marking



Fuse rating



Protective earth



Dangerous voltage / risk of electrocution



To be disposed off separately from household waste



Equipotential terminal



Refer to the instruction manual / brochure



General safety sign



Housing protection classification



Catalog reference



serial number



For indoor use only

Remark

the equipotential terminal allows the connection of an equipotential conductor to bring all the different elements of an electro-medical system to the same potential.



Do not sit on the back rest



Do not sit on the leg rest



You must not use the medical device for transfers



Do not sit on the armrests



Warning

You must not remove the pictograms and warning signs provided by the manufacturer! The manufacturer disclaims all responsibility in case of removal of these signs.

Applied parts

The applied parts according to standard EN 60601-1 are:

- PVC upholstery
- Leg rest
- Gynecological examination stirrups (accessories)

Electrical power supply



The equipotential terminal must not be used as a protective earth connection under any circumstances.



Equipped with the optional electrical outlets, 240-01xx, this medical device is electrical class 1. It must be connected to a grounded power supply system only.

Electromagnetic interference

This Promotal medical device was designed and built to minimize electromagnetic interference with other equipment. If interference is, however, observed, you must remove the apparatus causing the interference from the room and/or plug it into an isolated circuit.

Unpacking precautions

Medical device delivered on a wooden pallet

The medical device positioned on a wooden pallet may be easily moved using a forklift truck, as long as this is used correctly. Before transportation,

ensure that the forklift truck is correctly positioned in relation to the pallet, and that the unit is stable.

Storage conditions

Room temperature:

-15 °C to +60 °C (+5 °F to 140 °F)

Relative humidity :

10 % to 90 % (without condensation)

Atmospheric pressure :

0.5 bar to 1.05 bar (500 hPa to 1050 hPa)

All storage must be carried out in accordance with the following recommendations:

- Clean, aired and temperate area.
- Medical device stored in an area sheltered from bad weather and direct sunlight.
- Dry room.
- Medical device protected from shocks.

- Do not store in an area subject to frequent passage.

- Do not stack material.

- Keep in its original packaging until the final destination.

Conditions of use

- Dry and temperate area.
- Maximum altitude : 2000 m

- Temperature 10 to 40° C

- Relative humidity 75% maxi.

Unpacking and Installation


Step by step

1) During unpacking, remove all staples and remove the cardboard packaging carefully.

Caution : be careful with cutting tools, as fragile parts of the medical device (*covering, plastic housing, etc.*) may be near them.


2) If possible, transport the medical device on its pallet up to the final place of use.

The medical device is adjusted to a resistant position for transportation.

see: as indicated in the diagram 

3) Cut the two green bands that fix the chair to the pallet..

4) Next, take the chair off the pallet.



Caution

Four people are required to remove the medical device from the pallet.

Check

Having unpacked the medical device, follow these steps:

- 1) Check the delivery documents to ensure that the delivery is complete.
- 2) Check the external components for any

damage during transportation.

- 3) Check that the packaging contains the medical device, accessories and Options, the supply cable (*if electrical MD*) and the User Guide.

Remark:

Authorised EC Representative

Within the European Union, all problems, complaints or questions should be addressed to:

Promotal
22, rue de Saint-Denis de Gastines
53500 Ernée,
FRANCE
Telephone : + 33 (0)2 430 517 76
Fax : + 33 (0)2 430 572 00

Cleaning protocol

Warning

It is vital to read the user's manual thoroughly before manipulating this Medical Device. The equipment should only be used for its intended purpose as described in our documentation. Installation and connection must only be carried out by qualified personnel. The electrical components (cylinder, box,

control handle, battery, adapter, etc.) must not be opened under any circumstances. PROMOTAL shall not be held liable for any damage resulting from non-compliance with these instructions.



Any modification to the medical device without written authorization from the manufacturer is forbidden.



Caution

Only accessories designed and provided by Promotal for this medical device are authorised for use.



Caution

This medical device is not intended to be cleaned in a washing tunnel.

Cleaning/Disinfecting

This medical device must be regularly cleaned using the appropriate detergent products and disinfected using bactericidal, virucidal and fungicidal disinfectants.

A mild detergent such as soapy water can be used for routine cleaning of upholstery, stainless steel, aluminium or painted surfaces, plastic parts and control components, followed by effective rinsing and thorough drying. Detergents and disinfectants designed for use with medical appliances, such as those containing quaternary ammonium compounds, hydrogen peroxide, ethanol, chlorine compounds, etc. can be used on our medical devices provided that:

- The concentration prescribed by the suppliers of such products are complied with;
- The application conditions (contact time, quantity used, temperature, rinsing, etc.) are complied with;
- The supplier's instructions state that the detergent-disinfectant used is suitable for use with:
 - PVC, ABS, Polyamide, Polyurethane, Polypropylene
 - Epoxy-coated metal surfaces
 - Stainless steel or aluminium metallic surfaces.

Warning:

- Solvents are strictly prohibited.
- The use of abrasive powders or any other abrasive product should be avoided.
- High-pressure cleaning is forbidden.



Under no circumstances shall Promotal be held liable under warranty for any damage caused by non-compliance with the use instructions for a detergent-disinfectant.



Intended purpose

This medical device is designed for use in professional premises only:

- Medical office
- Health establishment

This medical device must not be installed in domestic premises.

This medical device is intended for general or specialist medical examinations and acts. Only health professionals (doctors, caregivers) and/or qualified members of the technical staff can handle or use the device.

Patients and companions must not intervene directly on the device. They must not have access to the device's controls.


During the consultation, supervision of children is the responsibility of the parents or of the practitioner.



This medical device must never be used as a healthcare bed.

This medical device must not be used in the presence of anaesthetics flammable in contact with air or nitrous oxide to avoid all risks of explosions.

Patient weight capacity



Protection against penetration of liquids


- IP X1

Characteristics


- Height 84 cm
- Steel frame with white epoxy finish
- Headrest adjustable by means of a gas lift
- Manual seat tilt (7°)
- Removable leg rest
- Upholstery Width 72 cm
- Seamless removable upholstery
- Upholstery with M1 fireproof coating
- Paper roll holder,
- Maximum length 50 cm – ø 20 cm (head end)
- Pull-out stirrups
- Removable plastic tray




Dimensions / Installation precautions



235-10 / 235-15



235-10 / 235-15 + Option 272-01




The table is delivered dismantled.
The packaging contains 2 boxes :
the frame


the upholstery, the paper roll holder.






Laying out the sleeping surface
One piece upholstery





 Ensures good clamping with each use

Legrest upholstery



**Caution**

When this medical device is equipped with option 240-01xx, the power supply cable may represent a fall risk. Be aware of the cable when walking around the device once it is connected to the mains.

**Safety note**

Unplug the power cord before moving the medical device.

**Warning**

Connecting equipment to the multiple socket outlet effectively leads to creating a Medical Electrical System and the result can be a reduced level of safety.

**Safety note**

Table must be positioned in normal use such that the power cords can be easily accessed to unplug.

**Warning**

No modification of this equipment is allowed.

**Caution**

To completely isolate the table from electrical mains supply, both power cords must be unplugged.

When using high frequency surgical devices or endocardial catheters:


- Use non-conductive material to insulate patient from metal portions of the table.
- Consult operating instructions for the device before using in conjunction with drawer or upholstery heater.

Failure to comply may result in electric shock or burns to patient.

Using the Kompact

Backrest tilt adjustment

With one hand, operate the lever (L1), and with the other hand, lower or raise the backrest according to the desired inclination.

 The patient should not lean on the backrest during adjustment.





Trendelenburg position


Principle

A. Lift the seat until it clicks into place.

B. With one hand hold the seat and with the other bring the lever (L) backwards. Replace the seat on the base.



 **Caution**
Make sure the device is locked before the patient gets on the table.

 **Caution**
of pinching when lowering the seat.

Leg rest

Using the leg rest



Caution

Maximum permissible load : 25 kg



Footstep


Using the footstep



Using the Kompact


Paper Roll

Paper roll installation



Using the paper roll

Cover the upholstery with paper before use.



Caution follow the indications provided well obtain a greater longevity of the material.

Gynaecology debris tray

The tray is removable for easy cleaning, and can be pulled out or pushed back under the seat.


In order to prevent the tray from falling out, the table is equipped with a retractable stop.



Using the gynaecology debris tray

A: Pull out

Push back on the retractable stop. Hold it down while pulling out the tray:



B: Push back in

Simply push the tray back into place.





Using the Kompact

Integrated stirrups

Installing the integrated stirrups

- 1) Pull the stirrup out, then unfold.
- 2) Lift the stirrup slightly, then move it left or right as desired.
- 3) Release stirrup to lock in position.



Be sure the stirrups are locked in place before using. The stirrups will not support the patient's entire weight.

The levellers

Adjustment of the levellers

Only concerns the 2 rear levellers.


Screw or unscrew the leveller (P) to the required height.





Drawers

The table is equipped with:

- 2 Pass through side drawers (accessible from both sides)
- 2 Front drawers with soft/self-close system



 **Caution**
Maximum permissible load : 8 kg

 **Caution**
Maximum permissible load : 5 kg

Using the Kompact

Braking


(only on models 235-15)

Using the brakes

2 pedal positions

1 : Press on the green foot pedal : All 2 castors are free.

2 : Press on the red foot pedal (brake position) : all 2 castors are locked.



Caution
Ensure brakes are fully engaged
before installing a patient.

Rear tray for storing paper rolls *(239-10)*

Set up


(Follow the instructions in the installation guide provided with this accessory).




2 removable plastic trays / Side drawers *(239-15)*

Set up

Open the drawer all the way and insert the plastic tray.



 Only accessories designed and provided by Promotal for this medical device are authorised for use.

Pair of stainless steel rails 25x10 (278-01)

VXM59820 (x4)



Remove plastic end caps



install the pair of rails


(Follow the instructions in the installation guide provided with this accessory).



Sliding clamp for 25 x 10 rail (reference 878-12)

Installing the clamp


878-12




Using the clamp




Check that the parts are secure before use.




 Only accessories designed and provided by Promotal for this medical device are authorised for use.

Rotating clamp for 25 x 10 rail (reference 879-10)



Installing the clamps

Insert the clamp on the rail.



Tighten the screws before use.

2 hook IV pole


(reference 985-01)

Installing the IV pole




Caution
⚠ Clamps are necessary to install these accessories.
(réf. : 878-12)

⚠ Ensures good clamping with each use



Auto-blocking I.V. pole

(reference 2985-01)



Installing the IV pole

Unscrew the screw handle.

Insert the I.V. pole and rescrew.


Caution
⚠ Clamps are necessary to install these accessories.
(réf. : 878-12)

⚠ Ensures good clamping with each use

Adjusting the height of the IV pole

The IV pole has 1 sliding stem.
Use the screw to adjust the height of the lower stem
and the bolt to adjust the upper stem.



 Only accessories designed and provided by Promotal for this medical device are authorised for use.


Arm Rest


(reference 236-01)

Adjusting the arm rest


The standard arm rests can be adjusted in height (H) and can be positioned in all directions due to the ball

and socket joint underneath the armrest. To adjust the armrests, loosen the screws (VS) and position as desired.



 Ensures good clamping with each use





Caution
 Clamps are necessary to install these accessories.
(réf. : 878-12)


Pair of legrests GOEPEL

(reference 245-01)

Installing the legrests






 Ensures good clamping with each use


Caution
 Clamps are necessary to install these accessories.
(réf. : 878-12)


Lamp holder (reference 295-01)

Fixing of the lamp on the support



 Ensures good clamping with each use




 Ensures good clamping with each use


Options

Adjustable paper-roll holder (reference 2056-02)


Installing the paper roll holder



1. Remove the standard paper roll holder.



2. Replace the paper roll holder.




Installation of the paper roll

The paper roll is positioned on a mobile axis with adjustable hooks at each end.

1. Unscrew the screw handle (VS).
2. To remove the paper roll holder 1/2 slightly (PR).
3. Insert the paper roll.

Cover the upholstery with paper before use.



Leg rest flush with seat


(reference 272-01)

Using the leg rest



Caution

The leg rest should only be used to support the patient's legs. Under no circumstances should it support the person's entire weight. Accordingly, the leg rest must be pushed in when the patient is getting on or off the table. We decline all liability in the event of non-compliance with this clause.



Caution

Maximum permissible load : 25 kg



Caution

Ensure that the leg rest is correctly positioned before use.

Additional mains plug and equipotential terminal (reference 240-01EU / 240-01UK / 240-01US)

Available only in EU, UK & US standard



Safety note

By connecting a medical device to the additional mains plug, the unit becomes an Electro-medical system according to the standard EN 60601-1. The user must ensure that the EM system is in conformity with the standard EN 60601-1 (article 16).



Caution

the supply cable for this medical device may represent an obstacle and cause falls. Do not forget its presence when moving around the device once connected to the mains supply.



Safety note

Unplug the power supply cable before moving the MD.


Using the additional mains plug

- Characteristics : 120 V or 230 V depending on the country
2,5 A maxi



The power consumed in a single event must not exceed 600VA.

Options (supplemental power socket) consists of 2 sockets, one on each side of the base of the couch.



The electrical sockets are equipped with hinged covers.

When these covers are closed, the examination couch is protected against sprays of water (IP X4 according to European standard 60 529). If one of the covers is open (i.e. each time one the sockets is in use), this protection is no longer present.

Equipotential terminal



This medical device is of electrical class 1. It must only be connected to a supply network equipped with protective earth.



Remark

the equipotential terminal allows the connection of an equipotential conductor to bring all the different elements of an electro-medical system to the same potential.



The equipotential bonding conductor must never be used as the earth connection for protection.

Lifespan of the Medical Device

This device is designed for a 10 year lifetime (except wear parts) under normal use conditions. This lifetime may vary according to frequency of use.

CHECK THE GENERAL CONDITION OF THE DEVICE AT LEAST ONCE A YEAR.

The wear parts are:

- The upholstery.
- The gas spring

Promotal recommends replacing the wear parts after 5 years use maximum.

For all interventions, contact your usual dealer, indicating the Serial number of the device.



Material warning

Within the framework of maintenance, only the installation of components designed and provided for this MD by Promotal is authorised.

Compulsory / specific maintenance

Once a year, ensure that the following checks are carried out by a qualified technician (contact your dealer):

- Ensure that all screws are correctly tightened.
- Check the fixings of articulated parts.
- Ensure that the structure has not been deformed.
- Ensure the fuses are in good condition.
- Check that the power cable has not been cut or damaged.
- Check the different connections (excessive play, noise...)
- Note this information and the control date in the maintenance record.

Once a month :

- Check the correct operation of the safety devices (cf page 16).

Medical device end of service life

Your dealer is responsible for the recovery and end of life treatment of this device.

If necessary, do not hesitate to contact Promotal. We can propose solutions to treat this equipment in the best conditions.



Material warning

Upholstery whose coating is torn no longer provides an effective anti-bacterial barrier and must be replaced without delay.



Material warning

DO NOT dismantle THE DEVICE

If a fault is detected, immediately contact your dealer or the dealer's technical department (cf. Warranty Chapter) for a complete diagnosis. If you have a doubt, do not use the device.

Fuse replacement

(Only if option 240-01EU or 240-01UL or 240-01US is present)


1. Remove all power to unit by unplugging unit's power cord.
2. Using a flat tipped screwdriver, gently release the plastic spring clip (A) on both sides of fuse drawer (1) then pull fuse drawer from fuse housing (2).
3. Pull both fuses (3) out of fuse drawer (1) and inspect. Check the fuses for any indication that they have blown ; i.e. burnt look, fuse cord melted through, etc. Discard fuses (3) if blown and replace them.



EQUIPMENT ALERT

Use fuses of the same voltage rating, amperage rating, and type. Failure to do so could result in damage to the equipment.

4. If necessary, obtain new fuse(s) (3). The replacement fuse(s) must be a 250 VAC, 6.3 amp, IEC 127 rated, 5 x 20 mm, Type T "Slo-Blo".
5. Insert fuse(s) (3) into fuse drawer (1).
6. Insert fuse drawer (1) into fuse housing (2) until fuse drawer snaps into place (both side of use drawer are locked into fuse housing).
7. Plug in power cord to equipment. If fuse blows again, call your Promotal distributor.



Warranty Information

Warranty

Promotal warrants, to the original purchaser, products manufactured by Promotal and components to be free from defects in materials and workmanship for a period of two (2) years¹ from the date of purchase.

Obligations

Promotal will, at its discretion and expense, replace defective parts reported to Promotal within the applicable warranty period, and which, upon examination by Promotal, prove to be defective.

In accordance with CE regulation, Promotal's distributor is responsible for after sales service during and after the warranty period.

Exclusions

- This warranty does not extend to :

- (1) Spare parts and consumables.
- (2) Travel and labour expenses.
- (3) Breakdowns due to improper use, manifest neglect, or to moving the device.
- (4) Equipment whose original characteristics have been changed by the user.
- (5) Control units, hydraulic and electric jacks if opened by the user (seals broken).
- (6) Damage, breakdowns, failures, or defects attributable to outside causes (lightning, electric surges, floods, natural catastrophes, impacts, etc.) or to the presence of foreign objects.
- (7) Damage caused by improper hook-up or by the power supply, damage caused by corrosion or by gradual deterioration of the product.
- (8) Indirect damage related to loss of use and penalties generated by poor performance.
- (9) Aesthetic damage incurred by the outside parts of the equipment, which does not hinder proper operation, such as scratches, chips, and scrapes.
- (10) Devices whose serial number has been rendered illegible or has been changed or removed

Exclusive Remedy

Promotal's only obligation under this warranty is the replacement of defective parts. Promotal shall not be liable for any direct, indirect, special, incidental, exemplary, or consequential damages or delay, including, but not limited to, damages for loss of profits or loss of use.

No Authorization

No person or firm is authorized to create for Promotal any other obligation or liability in connection with the products.

THIS WARRANTY IS PROMOTAL'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED.

¹ Upholstery sets are warranted 1 year from manufacturing defaults.



CE DECLARATION OF CONFORMITY

We,
Promotal

22, rue de Saint-Denis de Gastines
53500 Ernée
FRANCE

Hereby declare in sole responsibility that the medical device :

Model number : **Kompact 235-10 / 235-15**

Description : **Examination table**

is a medical device, class I,

we certify that the above-mentioned product meets the following requirements :

directive 93/42/EEC of the Council of the 14/06/93

(Medical device / appendices I and VII / decree 95-292 of the 16th March 1995),

all European standards applicable to it,

and herefore carries CE marking.

Year in which CE marking was first affixed (*Medical Device*) : **2020**

Ernée,
1 july 2020

A handwritten signature in black ink, appearing to read "R. Mouradian", with a long horizontal flourish extending to the right.

Rudolf MOURADIAN
President



CE

PROMOTAL - FRANCE
www.promotal.com
DIC235-10_2720EN

