

GYNAECOLOGICAL TABLE







THE EFFICIENT AND EASY TO USE CHAIR



COMFORT

- > 1 actuator
- Gas-spring assisted backrest Bilateral receptacle for optional colposcope holder Combined motion of the height and backrest
- > Foot control



ACCESSIBILITY

- Electric variable height from 46 to 106 cmCutout in the base for more accessibility to patient
- > Versatility of the table
- (gynecology and general examination with the cushion as accessory)
- > 3 models to adapt to various practices (stirrup, goeppels, footrests)



SAFETY

- > Capacity 170kg
- > Emergency position 11° > Equipotential bonding conductor: enables to equalise the potentials between dicerent electromedical devices





HYGIENE

- > Removable formed seamless upholstery Steel frame structure with white epoxy finish
- Paper roll holder
- Sliding and removable ABS tray on rails
 Waterproof seatrest
- > Casters with brakes (option)











STANDARD FEATURES

Electric height adjustment from 46 to 106 cm

Combined height and backrest movements (1 actuator)

Gas-spring assisted backrest adjustment

Equipotential terminal

Foot control

Adjustable levellers under the base

Sliding and removable plastic pan on rails (6 litre capacity)

Bilateral anchor points for the colposcope holder 's accessories

Steel frame structure with white epoxy finish

Removable formed seamless upholstery

Paper roll holder (max. length 50 cm - ø 12,5 cm)

170 kg maximum patient weight capacity (tested at a dynamic weight of 210 kg and a static weight of 680 kg)

11610-01 Stirrups

11614-01 Goeppels

11615-01 Footrest

OPTIONS*

Plant installed, must be ordered with initial order

ACCESSORIES




NET WEIGHT: 95 KG GROSS WEIGHT: 135 KG PACKAGE VOLUME: 0.92 M³


CHOICE OF M1ANTIBACTERIAL FIRE RETARDANT UPHOLSTERY






^{*} The Elegance range is subject to the application of a price supplement. Please refer to the price list.













WWW.PROMOTAL.NZ WWW.PROMOTAL.COM.AU