







A couch that is both easy to store and to install





STANDARD FEATURES

Epoxy polyester coated steel frame

Interlocking feet

Adjustable levellers

Flat upholstery, M1 fire retardant covering with antibacterial protection (depending on the color, see colour chart)

Adjustable backrest with trammel

Patient weight capacity 135kg

HYGIENE

- All the outer parts are accessible for cleaning
- Covering with antibacterial protection (depending on the color, see colour chart)

ACCESSIBILITY

• Height with upholstery 75 cm

SAFETY

- M1 fire retardant upholstery (depending on the color, see color chart)
- 2 adjustable levelers for an optimal stability
- Fully welded frame
- Capacity 135 kg





Net weight 32 kg Gross weight 34 kg Volume 0.210 m³





