



ICONA

gMotio



Updated: February 2022

PROMOTAL

gMotio

The gMotio gynaecology couch was designed following 2 years of study with health professionals in Europe and the Middle East. It has been tested by around 20 doctors and 100 patients. gMotio has been designed to provide effective and comfortable treatment, both for the practitioner and for their patients. With its extensive experience, Promotal has created a synergy between two featured products: the Elite gynaecology couch and the eMotio examination couch. gMotio is already considered a must-have gynaecology couch.





EQUALITY OF TREATMENT

MAKE LIFE EASIER FOR YOUR PATIENTS WITH REDUCED MOBILITY WITH:

- > Starts at 57cm from the ground
- > Electrical support accompaniment of all patient movements
- > Electrical moving footrests (depending on model)



EFFICIENCY

IMPROVE THE FLOW OF YOUR EXAMINATION WITH:

- > An integrated leg rest extension for a trendelenburg position (opt)
- > 3 programmable positions (9 as accessories)
- > Autonomous movements: Press once to change positions
- > Intuitive control and its already configured patient start position
- > Integrated lighting in the table (depending on model)
- > Standard seat rails



ERGONOMICS

ALWAYS PROVIDING THE BEST POSITION FOR YOU AND YOUR PATIENT:

- > Precise work on the seat design and foam density
- > Variable height from 57 to 118,5cm (with tilt)
- > Precise and adjustable tilt, easy speculum placement
- > Controls at the right place on the backrest and on the feet
- > Built in retractable wheels for short trips



SAFETY

EASY CLEANING OF YOUR COUCH:

- > Seamless upholstery
- > Removable seat and backrest for easy cleaning
- > Cataphoresis treatment of the structure for longevity and hygiene of the product
- > Removable ABS tray
- > Electrically activated wheels by sustained contact

SAFEGUARD YOUR EXAMINATIONS:

- > Safety sensors
- > Safety battery
- > Dual control

gMotio



ELECTRICAL FOOT RESTS + INTEGRATED LIGHTING

— Ref 15640-01 (pictured above)

ELECTRICAL FOOT RESTS

— Ref 15630-01

Model presented with accessories



STIRRUPS ON CLAMPS

— Ref 15620-01



LEG RESTS ON CLAMPS

— Ref 15626-01



GOEPEL ON CLAMPS

— Ref 15624-01

— OPTIONS

— ACCESSORIES

gMotio

- OPTIONS
- ACCESSORIES



Integrated mains socket and equipotential terminal
(EUROPE, UK, US) Ref. 40400-1
(DK, IT, AUS, CH) Ref. 40400-2



Urology tray
(replacing the std tray)
Ref. 2061-30



Stainless steel tray
(replacing the std tray)
Ref. 15966-01



Leg rest extension
Ref. 15672-01



PVC Glass extension leg rest protection
Ref. 15907-01




Pair of coated stirrups
(add x2 clamps 878-12)
Ref. 12850 (+ 2 x 878-12)
Incompatible with 15630-01
and 15640-01



Pair of plastic leg rests
(add x2 clamps 878-12)
Ref. 840 (+ 2 x 878-12)
Incompatible with 15630-01
and 15640-01



Pair of covers for leg rests
Ref. 861



Pair of goepfels
(add x2 clamps 878-12)
Ref. 845 (+ 2 x 878-12)
Incompatible with 15630-01
and 15640-01



Single armrest
(add x1 clamp 878-12)
Ref. 2536
(+ 878-12)



Armrest cover
(for models 2536)
Ref. 2537



Footrest hoop cover
(for models 15630-01 & 15640-01)
Ref. 15861-01



Adjustable headrest
(iQuest 3055-20 and gMotio)
Ref. 11679-01



PVC Glass footrest protection
Ref. 15905-01



PVC Glass seat protection
Ref. 15906-01



Manual remote control
(9 programmable positions)
Ref. 15031-10



Self-locking serum holder rod 2 hooks
Ref. 2985-01 (+ 878-12)



Serum rod holder 2 hooks with stop
Ref. 985-01 (+ 878-12)



Splash guard
Ref. 12677-01



Colpo support. Zeiss d.38 for gMotio*
Ref. 15840-01

*suppresses autonomous movements



Colpo support. Pivoting Leisegang for gMotio*
Ref. 15810-01

*suppresses autonomous movements



Olympus right colpo support for gMotio*
Ref. 15850-01

*suppresses autonomous movements



Olympus left colpo support for gMotio*
Ref. 15851-01

*suppresses autonomous movements




Stainless steel backrest rail (pair) 25x10 mm
Ref. 15879-01



Sliding clamp for rail 25x10 mm, ø 16/18 mm
Ref. 878-12



Rotary clamp for rail 25x10 mm
Ref. 879-10



Lamp support for rail* 25x10 mm
Ref. 295-01

*to select a light, please refer to the specific documentation




Lamp
(Refer to the documentation)

gMotio

MODELS 15620-01/15624-01/15626-01 :
 NET WEIGHT : 177KG
 GROSS WEIGHT : 210KG
 PACKAGE VOLUME : 1.02M3

MODEL 15630-01 :
 NET WEIGHT : 181KG
 GROSS WEIGHT : 214KG
 PACKAGE VOLUME : 1.02M3

MODEL 15640-01 :
 NET WEIGHT : 182KG
 GROSS WEIGHT : 215KG
 PACKAGE VOLUME : 1.02M3



15620-01 / 15624-01 / 15626-01

15630-01 / 15640-01

UPHOLSTERY COLOURS

CA - TOMATO		CLASSIC RANGE
CB - FUCHSIA		
CC - PLUM		
CD - LEMON		
CE - PEACOCK		
CF - SAPHYR		
CG - MIDNIGHT BLUE		
CH - GLACIER		
CI - DARK		
57 - ANTISTATIC		
CJ - GREY		ELEGANT RANGE*
CK - PLATINUM		
EA - DEEP GREEN		
EB - VEGETAL		
EC - STORM		
ED - SILVER BLUE		
EE - PAPRIKA		
EF - TWILIGHT		
EG - CAMEL		
EH - CORK		

* The Elegance range is subject to the application of a price supplement. Please refer to the price list.