



STANDARD • HOSPITAL • LOUNGE





In 10 years, the Quest examination couch has satisfied more than 15,000 healthcare professionals. General practitioners, cardiologists, dermatologists and physiotherapists from over 70 countries have come to appreciate

the French quality of this robust, comfortable and easy to access examination couch. Five years ago, Promotal conducted the largest survey ever performed in Europe and the Middle East. This 3-year study, carried out in collaboration with more than 200 health professionals, revealed 3 basic needs: equal treatment of people with mobility limitations, efficiency in carrying out medical examinations and user-friendly equipment. iQuest is Promotal's new revolution. We have designed an examining couch with a sleek, lightweight design that will improve your comfort and that of your patients thanks to the low position that is the lowest on the market, the programmable positions and the autonomous movements that are 2 times faster.















EXAMINATION ACCESS

HIGH **POSITION**

TRENDELENBURG

REVERSE TRENDELENBURG



EQUAL TREATMENT

MAKE THINGS EASIER FOR YOUR PATIENTS WITH MOBILITY LIMITATIONS:

- **>** Low 46 cm position provides easier access for patients
- > Patients can use the raised backrest for support when getting on the couch
- > Safe, smooth electrical movements

EXAMINE ALL YOUR PATIENTS EASILY:

- > Accommodates up to 250 Kg
- > Superior comfort with the Lounge upholstery



EFFICIENT DESIGN

YOUR EXAMINATIONS WILL GO SMOOTHLY THANKS TO

- > Movements that are 2 times faster
- > The intuitive position of the patient on the Lounge upholstery
- **>** The autonomous movements: Press once to change positions
- > The intuitive controls and the pre-programmed patient access position



USER-FRIENDLY DESIGN

CHANGE POSITIONS EFFORTLESSLY:

- > All-electric examination couch: Height, backrest, Trendelenburg / reverse Trendelenburg
- > Variable height from 46 to 93 cm
- **>** 2 programmable positions, 100% configurable
- > Intuitive controls built into the backrest (Lounge model)

IMPROVE YOUR PATIENTS' COMFORT THANKS TO:

> The superior comfort Lounge upholstery



3 MODELS, 3 NEEDS.

iQuest

STANDARD

Ref. 3050-40 With retractable wheels Ref. 3050-50 Upholstery width 65 cm Ref. 600150-65 **Upholstery width 75cm** Ref. 600150-75

iQuest

LOUNGE

Ref. 3055-20










Face cut-out 65cm **Ref. 2231-65** 75cm **Ref. 2231-75**

upholstery Ref. 600000-01

Battery Ref. 3092-01



Circuit breaker Ref. 3058-01



Ref. 3031-01



Integrated plug Equipotential EU-US-UK Ref. 30400-01 Ref. 2059-01





Clamp Ref. 166



Pair of clamps on rail

Ref. 878-12

Pair of legrests Ref. 840

Armrest

(one) **Ref. 2536**



Pair of legrest covers Ref. 861



Pair of stirrups Ref. 850



IV pole with 2 hooks Ref. 985-01



Auto-locking IV pole with 2 hooks Ref. 2985-01



Armrest cover (one) **Ref. 2537**



Lounge headrest Ref. 30679-01



Flat upholstery headrest Ref. 1102-10



Adjustable paper roll holder Ref. 2056-02



Backrest extension Ref. 2055-01



Wired foot control Ref. 3051-10



Plastic cover (foot end) for flat upholstery
65cm Ref. 600150-65C 75cm Ref. 600150-75C



Plastic cover Plastic cover (head end) (foot end) Ref. 3056-01 Ref. 3056-02



Pair of head rails


Ref. 30878-XXT

Pair of foot rails Ref. 30878-XXP



Ref. 3050-60

iQuest



COLORS OF UPHOLSTERIES & COUNTERPLATES



*Colors benefiting from a matching counterplate. The other colors are equipped with a white counterplate





