



Quest

STANDARD • HOSPITAL • LOUNGE



iQuest

In 10 years, the Quest examination couch has satisfied more than 15,000 healthcare professionals. General practitioners, cardiologists, dermatologists and physiotherapists from over 70 countries have come to appreciate

the French quality of this robust, comfortable and easy to access examination couch. Five years ago, Promotal conducted the largest survey ever performed in Europe and the Middle East. This 3-year study, carried out in collaboration with more than 200 health professionals, revealed 3 basic needs: equal treatment of people with mobility limitations, efficiency in carrying out medical examinations and user-friendly equipment. iQuest is Promotal's new revolution. We have designed an examining couch with a sleek, lightweight design that will improve your comfort and that of your patients thanks to the low position that is the lowest on the market, the programmable positions and the autonomous movements that are 2 times faster.













ACCESS

EXAMINATION

HIGH POSITION TRENDELENBURG

REVERSE TRENDELENBURG



EQUAL TREATMENT

MAKE THINGS EASIER FOR YOUR PATIENTS WITH MOBILITY LIMITATIONS :

- > Low 46 cm position provides easier access for patients
- > Patients can use the raised backrest for support when getting on the couch
- > Safe, smooth electrical movements

EXAMINE ALL YOUR PATIENTS EASILY :

- > Accommodates up to 250 Kg
- > Superior comfort with the Lounge upholstery



EFFICIENT DESIGN

YOUR EXAMINATIONS WILL GO SMOOTHLY THANKS TO

- > Movements that are 2 times faster
- > The intuitive position of the patient on the Lounge upholstery
- > The autonomous movements: Press once to change positions
- > The intuitive controls and the pre-programmed patient access position



USER-FRIENDLY DESIGN

CHANGE POSITIONS EFFORTLESSLY :

- > All-electric examination couch: Height, backrest,
- > Trendelenburg / reverse Trendelenburg
- > Variable height from 46 to 93 cm
- > 2 programmable positions, 100% configurable

Intuitive controls built into the backrest (Lounge model)

IMPROVE YOUR PATIENTS' COMFORT THANKS TO :

> The superior comfort Lounge upholstery





Ref. 3050-40 With retractable wheels Ref. 3050-50 Upholstery width 65 cm Ref. 600150-65 Upholstery width 75cm Ref. 600150-75

iQuest LOUNGE

Ref. 3055-20 With retractable wheels





MODEL 3050-40 (UPHOLSTERY NOT INCLUDED) NET WEIGHT : 100KG GROSS WEIGHT : 130KG PACKAGE VOLUME : 1,14M³

MODELS 3050-40/3050-60 (UPHOLSTERY NOT INCLUDED) NET WEIGHT : 105KG GROSS WEIGHT : 135KG PACKAGE VOLUME : 1,14M³

> UPHOLSTERIES 600150-65/600150-75 : NET WEIGHT : 21KG GROSS WEIGHT : 25KG PACKAGE VOLUME : 0,27M³

MODEL 3055-20 (UPHOLSTERY INCLUDED) NET WEIGHT : 127KG GROSS WEIGHT : 160KG PACKAGE VOLUME : 1,64M³


COLORS OF UPHOLSTERIES & COUNTERPLATES





* The Elegance range is subject to the application of a price supplement. Please refer to the price list.

**Colors benefiting from a matching counterplate. The other colors are equipped with a white counterplate










HOSPITAI

LOUNGE



PROMOTAI

22 rue de Saint-Denis-de-Gastines - **53500** Ernée, FRANCE Phone : + 33 (0)2 **43 05 68 95**



68

460 - 930

68

460 - 930

15

ç

15

174

N⊙M≜D⊒° www.agence-nomade.re